



## Vascularperspectives

innovations for the interventionalist

**BEO 1.5mm RBP 24atm**

NC OTW multi-tasking balloon

**BEO 1.5 NC OTW Balloon:**

- Bi-flex tip combining two essential elements for crossing heavily calcified lesions and tortuous vessels or for crossing stent struts into side-branches:
  - 1) Distal red 2mm section is highly flexible and provides superior over the wire trackability and lesion entry
  - 2) Yellow 3mm section provides stiff-tip-support to push through tight lesions or stented side-branches.
- High RBP with low compliance (< 12% at 24atm)
- 160cm translucent shaft enabling wire visualization during catheter exchange
- Long enough for retrograde approaches to CTO's

**BEO 1.5 NC OTW Balloon Features**

- 1 5mm Bi-flex tip
- 2 Soft tip (red=2mm)
- 3 Stiff support (yellow=3mm)
- 4 RBP 24atm, 1.5mm at 8atm and 1.68mm at 24atm
- 5 160cm translucent shaft

**Ordering information****OTW BEO 1.5mm NC High Pressure Balloon**

Cat. Number	Balloon Diameter	Length	RBP	Min. Guiding
150-010-034	1.5mm	10mm	24	5F
150-015-034	1.5mm	15mm	24	5F
150-020-034	1.5mm	20mm	24	5F

# Vascularperspectives

innovations for the interventionalist

Parkview House  
127 Styal Road  
Heald Green, Cheadle  
Cheshire SK8 3TG, UK

Tel: 0161 491 0585  
Fax: 0161 282 6540  
Web: [vascularperspectives.com](http://vascularperspectives.com)  
Email: [info@vpmed.co.uk](mailto:info@vpmed.co.uk)